



## BARIATRIC

REORDER NO. 10233 – 22"×18"  
REORDER NO. 10234 – 24"×18"

## BARIATRIC HD

REORDER NO. 10290 – 26"×20"  
REORDER NO. 10291 – 28"×20"  
REORDER NO. 10292 – 30"×20"  
REORDER NO. 10293 – 32"×20"

## BARIATRIC WHEELCHAIRS

### WITH ELEVATING LEG RESTS AND DETACHABLE DESK ARMS

- Powder-coated steel frame with dual cross brace and detachable desk-length armrests
- Black, padded, vinyl (10233-10234) or nylon (10290-10293) upholstery
- Dual-axle provides adjustment from standard height to hemi level (10233-10234)
- Tool-free, adjustable, elevating leg rest features extra-wide foot plate and padded calf
- 24" composite real wheels with solid rubber tires
- Push-to-lock brakes
- Includes built-in chart pocket
- Limited Lifetime Warranty
- HCPCS: K0007, E0990 OR K0195

## WEIGHT CAPACITY

- **500 lb.** (10233-10234)
- **600 lb.** (10290-10292)
- **700 lb.** (10293)



**BARIATRIC**  
10233 - 10234



**BARIATRIC HD**  
10290 - 10293