

INVACARE®
9000 TOPAZ™
WHEELCHAIR



The **Invacare® 9000 Topaz™ Wheelchair** is designed specifically to suit the unique needs of bariatric patients. The Topaz wheelchair combines the strength and durability of a carbon steel frame with the comfort of heavy-duty front riggings and arms for proper positioning. With a variety of seat widths and depths, frame reinforcements that include gusseted sides, caster journals and dual heavy-duty crossbraces, this chair provides not only ease in maneuverability but proven stability.

Features

Available in vinyl or nylon upholstery

Durable low-maintenance, carbon steel frame is long-lasting

Reinforced frame includes gusseted sides, caster journals and dual heavy-duty crossbraces

Multiple frame color variations available for customizing

Heavy-duty front riggings provide proper positioning for feet with footrest and elevating legrest

Heavy-duty arms designed to accommodate weight capacity and provide proper positioning

Dual axle positions allow for quick and simple seat-to-floor height variations, 17.5" and 19.5"

Heavy-duty package option for 1,000 lb. weight capacity

Heavy-duty triple inner liners keep seat and back from stretching

HCPCS Code: K0007



1. Dual axle positions
2. Heavy-duty front riggings
3. Heavy-duty arms
4. Bariatric sizes
5. Seating available
6. Builder options



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2014
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the
symbols ™ and ®.
All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 11-230
Rev. 04/14
131072

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
9TPZ	39"	32"	30"	16"	20"-30"(2" inc) x 18", 20"	15.5"-19.5"	16"-18"
	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
9TPZ	6", 8"	20", 24"	55 lb.	700 lb., 1,000 lb	97 lb.	35"x15"x 46.5"	3 Years

