



Yes, you can.®

INVACARE HME BED FRAMES OVERVIEW

Homecare and LTC Bed Systems

AGENDA



Yes, you can.®

 Category Overview

 Dedication to Patient Safety

 Supply Chain

 Product Overview



☰ CATEGORY OVERVIEW

INVACARE HOMECARE

- Split frame design
- Quick assembly
- Universal Bed Ends

INVACARE CONTINUING CARE

- Outstanding bed design
- Pioneer of low bed design
- Tool less assembly on CS Line
- Large style offering



DEDICATION TO PATIENT SAFETY



FDA Entrapment Guidelines

- ✓ The term "**bed rail entrapment**" describes an event in which a patient/resident is caught, trapped, or entangled in the space in or about the bed rail, mattress, or hospital bed frame.
- ✓ Invacare worked on the team that **defined the entrapment zones** and guidelines that the FDA publishes.

Entrapment within
the bed rail



Entrapment under
the bed rail



Entrapment between
the rail and mattress



Entrapment under
the bed rail at the ends
of the bed rail





DEDICATION TO PATIENT SAFETY



INVACARE HOMECARE

- ✓ Homecare beds are tested to [UL962](#)
 - Manufacturing Control
 - Reduced Shock Hazards
 - Platform Stability
 - Performance and Spec Validation
 - Lifecycle testing
- ✓ IVC Series carry an [IP66 Rating](#)
 - Water resistant connections to prevent fluid seepage
 - Increases infection control since it can be cleaned easier



INVACARE CONTINUING CARE

- ✓ DLX & CS series tested to [UL60601-2-38](#)
 - Includes entrapment testing
 - Eliminates Pinch Points
 - Increased scrutiny on electronics



SUPPLY CHAIN



HEMOCARE



IVC Series

- Stocked in OH, GA, CA
- Same day shipping



BAR Series

- Stocked in OH, GA, CA
- Same day shipping



Same day



CONTINUING CARE



DLX Series

- Built to order
- 3 Week lead time to ship
- Quick Ship Program available on <10 units



3 weeks



Quick Ship Program 72h (<10)



CS Beds

- Built to order
- 4 Week lead time to ship
- Quick Ship Program available (CA, GA)



4 weeks



Quick Ship Program same day
(CA - GA)



INVACARE BED PORTFOLIO



Yes, you can.®



HOMECARE BEDS



IVC Manual



IVC Semi Electric



IVC Full Electric



IVC Full Electric Low Height



BARIATRIC BEDS



BAR600 IVC



BAR750



LONG TERM CARE



820 DLX



SC900 DLX



Carroll CS3



Carroll CS5/CS7



Carroll CS9 FX600



Yes, you can.®

HME Beds





HME Bed Portfolio



Yes, you can.®



HEMOCARE BEDS

IVC Manual



Manual bed height adjustment



Manual body position adjustment

IVC Semi Electric



Manual bed height adjustment



Electric body position adjustment

IVC Full Electric



Electric bed height adjustment



Electric body position adjustment

IVC Full Electric Low Height



Electric bed height adjustment



Electric body position adjustment



Extra lower feature



INVACARE IVC Semi-Electric (5310IVC)



Meets UL962

Universal Bed Ends



Durable rubber pendant to prevent failures from dropping




Washable electronics



2 motor design to maximize rental cycles



Crank through motors in case of power failure





INVACARE IVC Full Electric (5410IVC)



Meets UL962

Universal Bed Ends



Durable rubber pendant to prevent failures from dropping



Washable electronics



3 motor design to maximize rental cycles



Crank through motors in case of power failure



INVACARE IVC Full Electric Low (5410LOW)



Universal Bed Ends



Meets UL962



Patented Assist Rail available to prevent pinch points and obstruction during Hi/Lo cycle




Low bed design to bring the patient lower to the floor



May help reduce injuries from falls¹



Electric Hi/Lo adjustment assists with transfers





HME Bed Portfolio



HEMOCARE BEDS (Electric)



	5310IVC	5410IVC	5410LOW
Description	Semi Electric	Full Electric	Full Electric Low Height
Motors	2	3	3
Height range	15" – 23"	15" – 23"	9.5" – 20"
Pendant	Rubber, 4-button	Rubber, 4-button	Rubber, 4-button
Bed Ends	Universal	Universal	Universal
Deck	Spring/Link Fabric	Spring/Link Fabric	Spring/Link Fabric
Heaviest Piece	65lbs	69lbs	69lbs
Warranty	2 years electronics 5 years on frame	2 years electronics 5 years on frame	2 years electronics 5 years on frame



HME Bed Portfolio



Yes, you can.®




HEMOCARE BEDS



	5310IVC	5410IVC	5410LOW
Description	Semi Electric	Full Electric	Full Electric Low Height
Head section	5000IVC	5000IVC	5000IVC
Bed Ends	5301IVC	5301IVC	5301LOW
Foot section	5890IVC	5490IVC	5490LOW

HME Bed Portfolio key points



🏠 HME Bed Portfolio key points



Double insulation



US-Based manufacturing



Thermal overload protection



Crank-through electronics



Motor Guard



Washable, color-coded electronics



Yes, you can.®

BARIATRIC Beds





Bariatric Bed Portfolio



Yes, you can.®



BARIATRIC BEDS

BAR600 IVC



BAR750





INVACARE BAR600 IVC Bariatric Bed



Meets UL962 and CSA certified

Split spring design with interchangeable bed ends for ease in set-up




42" by 80" sleep surface
(17% larger than conventional beds)



Bariatric 600lb total weight capacity



Easy release motor mounts





INVACARE BAR750 Bariatric Bed



Yes, you can.®

Bariatric 750lb total weight capacity

Split spring design for easy delivery and setup



Comes with bed ends, and half rails



Integrated battery backup



Expandable Length and Width

39" to 48" wide
80" to 88" in length



4 DC Actuators

2 hi/lo motors (one on head section/one on foot section)



Locking Casters

2 locking casters on HS
2 steering lock casters on FS





Bariatric Bed Portfolio



BARIATRIC BEDS



	BAR600 IVC	BAR750
Motors	3	4
Height range	14" – 22"	15.5" – 29.5"
Pendant	6-button	6-button
Bed Ends	Universal	Heavy-duty panels
Deck	Slat	Slat
Weight Capacity	600lbs	750lbs
Warranty	2 years electronics 5 years on frame	1 years electronics 5 years on frame



Yes, you can.®



THANK YOU