








Yes, you can.

INVACARE GROUP I MATTRESS

					
	Economy Foam	BARMATT	Solace Prevention	Softform Premier	Softform Glissando
HCPSC Code	E0271	E0272	E0184	E0184	E0184
Risk Level	Low Risk	Low Risk	Low - Medium Risk	Low - High Risk	Low - High Risk
Government Contract			VA National Standard #36C10G20D0018		

Mattress Construction

Material	High quality 312 coil innerspring design	High resilience foam designed for extended life while providing comfort and support	High quality specialized foam with three zones to help relieve pressure in the head, body, and heel sections	Unique contoured foam insert designed to provide excellent user comfort and optimize pressure reduction for users	Unique contoured foam insert designed to provide excellent user comfort and optimize pressure reduction for users
Cover	Durable 11-ounce vinyl cover, anti-bacterial, anti-static, acid resistant and waterproof	Durable 11-ounce vinyl cover, anti-bacterial, anti-static, acid resistant and waterproof	70 Denier Nylon cover and tri-laminate vinyl bottom cover for ease in cleaning	Four-way stretch fluid resistant polyurethane cover	Four-way stretch fluid resistant polyurethane cover
CetiPUR-US			✓	✓	✓

Mattress Features

Castellation Cuts		Solid foam block	The convoluted fortress cut design offers more surface area for resident support and relief	Premium HR foam cut into mushroom-shaped cells to allow superior patient immersion and pressure redistribution	Premium HR foam cut into mushroom-shaped cells to allow superior patient immersion and pressure redistribution
Therapeutic Zones	1	1	3	3	3
Density			1.5-1.8 lbs./ft ³	2.5- 3.7 lbs./ft ³	2.5- 3.7 lbs./ft ³
Anti-Shear and Friction Gliding Technology					✓
Firm Foam Pre-meter				✓	✓
Packaging			Compressed and Rolled	Compressed and Rolled	Compressed and Rolled

Item Numbers and Mattress Sizes WxLxH	5185 36"x80"x6" 5185XL 36"x84"x6"	BARMATT42 42"x80"x6.5" BARMATT48 48"x80"x6.5"	SPS1080 36"x80"x6" SPS1084 36"x84"x6" SPS2080B42 42"x80"x6"	IPM1080 36"x80"x6" IPM1084 36"x84"x6" IPM1080B42 42"x80"x6" IPM1080B48 48"x80"x6"	SRS2080 36"x80"x6" SRS2084 36"x84"x6" SRS2080W42 42"x80"x6"
Weight Capacity	350lbs.	750 lbs.	350 lbs.	500-650 lbs.	500 lbs.
Fire Code Regulations	16CFR1632/1633	16CFR1632/1633	16CFR1632/1633	16CFR1632/1633	16CFR1632/1633
Latex free	✓	✓	✓	✓	✓
Limited Warranty	1 Year	1 Year	1 Year Cover 5 Years Foam	4 Years Cover 8 Years Foam	4 Years Cover 8 Years Foam

INVACARE

GROUP II MATTRESS

	microAIR MA500	microAIR MA600	microAIR MA800	microAIR MA900
HCPCS Code	E0277	E0277	E0277	E0277
Risk Level	Low - Medium Risk	Medium - High Risk	Medium - High Risk	Medium - High Risk
Government Contract	VA Federal Supply Schedule #V797D-70085	VA Federal Supply Schedule #V797D-70085	VA Federal Supply Schedule #V797D-70085	VA Federal Supply Schedule #V797D-70085

Mattress Construction

Cover	Highly vapor permeable and oversized pliable quilted nylon cover	Highly vapor permeable and oversized pliable quilted nylon cover	Highly vapor permeable and oversized pliable quilted nylon cover	Highly vapor permeable and oversized pliable quilted nylon cover
Air Cells	18 air cells made from 70 Denier urethane coated nylon	18 air cells made from 70 Denier urethane coated nylon	18 air cells made from 70 Denier urethane coated nylon	18 air cells made from 70 Denier urethane coated nylon
Bolsters				Built in side bolsters
Cable Management	✓	✓	✓	✓
CPR Quick Release	✓	✓	✓	✓
Power Failsafe	Cell in Cell	Cell in Cell	2" foam base	2" foam base

Therapy Type

Alternating Pressure	2-in-1 Alternating Pressure Therapy with fixed cycle time of 10 mins	3-in-1 Alternating Pressure Therapy with adjustable cycle times ranging from 5, 10, 15, 20 mins	3-in-1 Alternating Pressure Therapy with adjustable cycle times from 3-95 mins	
Pulsation			Reduces the air flow every 60 seconds to 70% of the comfort level setting	
Low Air Loss	✓	✓	True Low Air Loss	True Low Air Loss
Lateral Rotation				20 and 40 degree lateral rotation with adjustable cycle times from 3-95 mins
Heel Slope	Natural Heel Slope	Heel Sense Technology	Heel Sense Technology	

Pump Features

Compressor	5 LPM	10 LPM		
Blower			1275 LPM	1275 LPM
Adjustable Comfort Level Settings	✓	✓	✓	✓
Auto Firm		✓	✓	✓
Audio and Visual Alarms	✓	✓	✓	✓
System lockout		✓	✓	✓

Item Numbers and Mattress Sizes WxLxH	MA500 36"x80"x8"	MA600 36"x80"x8"	MA800 36"x80"x8" MA800B42 42"x80"x8" MA800B4284 42"x84"x8" MA800B48 48"x80"x8"	MA900 36"x80"x8" MA900B42 42"x80"x8" MA900B48 48"x80"x8"
Weight Capacity	350 lbs	500 lbs	Up to 1,000 lbs	Up to 1,000 lbs
Latex free	✓	✓	✓	✓
Limited Warranty	1 Year Cover 1 Year Pump	1 Year Cover 1 Year Pump	1 Year Cover 1 Year Pump	1 Year Cover 1 Year Pump