

INVACARE®
TRACER® EX2
WHEELCHAIR



The **Invacare® Tracer® EX2 wheelchair** sets the benchmark in manual wheelchairs. A true dual axle position allows for repositioning the 24" rear wheels and 8" casters for adult and hemi seat-to-floor heights. Offset permanent arms offer full seat width up to 20" with a permanent arm, hemi-spaced footrest and full range of accessories. Incorporating the design and technology of the Invacare® 9000 series, the Tracer EX2 wheelchair has interchangeable components with the 9000 series chairs.

Features

- Durable, low-maintenance, carbon steel frame is long-lasting
- Dual axle positions allow variation of seat-to-floor heights
- Offset permanent arms available in 16", 18" and 20" widths
- 14-gauge crossbraces add strength and durability
- Urethane rear tires mounted on "no flex" wheels offer superior performance
- Low-maintenance, dual-sealed precision bearings
- Heavy-duty inner liners keep seat and back from stretching
- Upholstery is dual-embossed
- Meets RESNA standards for manual wheelchairs in product performance

HCPCS Code: K0001/K0002



1. Dual axle positions
2. Dual-embossed upholstery
3. Offset permanent arms
4. Seating options available



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 11-224
Rev. 03/13
130293



Model #	Overall Width (1" increments)	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Width	Seat Depth	Seat-to- Floor Height
TREX20	23.25"-28.25"	34"	31.88"	13"	20"	16"	17.5", 19.5"
TREX26	23.25"-28.25"	34"	31.88"	13"	16"	16"	17.5", 19.5"
TREX28	23.25"-28.25"	34"	31.88"	13"	18"	16"	17.5", 19.5"

	Back Height	Caster Height	Wheel Height	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
TREX20	16"	8"	24"	250 lb.	43 lb.	32"x13"x 38"	5 Years
TREX26	16"	8"	24"	250 lb.	43 lb.	32"x13"x 38"	5 Years
TREX28	16"	8"	24"	250 lb.	43 lb.	32"x13"x 38"	5 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.